



Testing of electrophoretic pumps in a brain tissue phantom

08/12/2021

Etienne Rognin (ecr43@cam.ac.uk)



UNIVERSITY OF
CAMBRIDGE
Department of Engineering





EPSRC IRC
Targeted Delivery
for Hard-to-Treat
Cancers

EPSRC

Engineering and Physical Sciences
Research Council

IfM | MANAGEMENT
TECHNOLOGY
POLICY

Overview of the testing rig



Absorbance images

- Absorbance A :


$$A = -\log_{10} \frac{\phi_1}{\phi_0}$$

- Beer-Lambert law:

$$A(\lambda) = \varepsilon(\lambda) \ell c$$





- ε Molar extinction coefficient

- Images: pixel-wise approach
 - First image of stack ($t=0$) as reference




Dynamic range



- Max and min concentrations depend on pixel depth.
- Other sources of limitation:
 - Sensor non-linearity
 - Sensor noise
 - Back-light non-uniformity
 - Scattering in the gel



1D diffusion of methylene blue (MB)



1. Blank gel
(reference)



2. MB ($3\mu\text{g/ml}$)
is poured on top




3. Absorbance recorded at
1 frame per minute




1D diffusion of methylene blue

After 1h




After 12 hours




1d diffusion: results

- Apply function fit to the full stack of images
- Function parameters give diffusivity coefficient
- For water, $D = 6.7 \times 10^{-6} \text{ cm}^2/\text{s}$.





3D diffusion (sharp needle)



Analysis after needle removal

- Axial symmetry assumed
- Projected absorbance of a Gaussian is also a Gaussian

$$A(y) = \varepsilon c_0 \sigma \sqrt{2\pi} e^{-\frac{y^2}{2\sigma^2}} \iff c(x, y) = c_0 e^{-\frac{x^2 + y^2}{2\sigma^2}}$$



Analysis over time

- In Fourier space, diffusion is a Gaussian filter:

$$\hat{c}(\mathbf{k}, t) = \hat{c}(\mathbf{k}, t = 0)e^{-k^2 Dt}$$


- The Fourier transform of a Gaussian is also a Gaussian, & vice versa:

$$c(x) = c_0 e^{-\frac{x^2}{2\sigma^2}} \Rightarrow \hat{c}(k) = c_0 \sqrt{2\pi\sigma^2} e^{-\frac{\sigma^2}{2} k^2}$$

- So, the scale of a diffused Gaussian evolves as:

$$\sigma(t)^2 = \sigma(0)^2 + 2Dt$$


- We can measure the diffusivity D




Diffusion from devices (retrodialysis)




1. Hot gel is poured



2. Wait 1 hour,
then connect to
syringe pump





3. Flush with 0.2 mL of MB



4. Start syringe
pump,
capture 1 frame/min



Diffusion from devices (retro-dialysis)



No flow

1 $\mu\text{L}/\text{min}$

waste

Concentration map

- De-noising
- Computed tomography assuming axial symmetry







Image extraction



Symmetry axis




Absorbance
After 700 min



Concentration
Note: inlet at 10
 $\mu\text{g/mL}$

Integrated mass


Sum of absorbance pixels gives the mass of MB in the cuvette (minus probe shadow).






Influence of flow rate

Péclet number in device:

$Pe = \text{Flow rate} /$
 $(\text{diffusivity} * \text{Length})$



Electrophoretic operation



Conclusion

Test rig can be used to validate probes in retrodialysis operation

- Accurate concentration fields of MB
- Different flow rates
- Other important parameters: pressure, connectors.
- This will be used to validate models and simulations and feed design iterations.

Current work

- Electrophoretic operation
- DC and AC modes, electrodes considerations.
- Testing of conductive elastomers in gels.

Acknowledgements

Chemistry

Dr Niamh Willis-Fox

Probe manufacturing

Tobias Naegele

Co-Investigator

Dr Ronan Daly

EPSRC

Engineering and Physical Sciences
Research Council



EPSRC IRC
Targeted Delivery
for Hard-to-Treat
Cancers



**UNIVERSITY OF
CAMBRIDGE**

Department of Engineering

IfM | MANAGEMENT
TECHNOLOGY
POLICY